










UCD DOPLM IMMUNOLOGY/FLOW CYTOMETRY ALTERNATE SPECIMEN TUBES

TEST	Approved tube types for laboratory testing
<ul style="list-style-type: none"> ▪ Protein Electrophoresis-Serum ▪ Immunofixation-Serum ▪ Anti-Mitochondrial Antibodies ▪ Anti-Parietal Cell Antibodies ▪ Anti-Smooth Muscle Antibodies 	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  SST </div> <div style="text-align: center;">  Red Top </div> </div>
<ul style="list-style-type: none"> ▪ RAST Allergy Testing <ul style="list-style-type: none"> ○ Specific IgE ○ Total IgE 	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  SST </div> <div style="text-align: center;">  Red Top </div> <div style="text-align: center;">  EDTA </div> <div style="text-align: center;">  Sodium Heparin </div> <div style="text-align: center;">  PST </div> </div>

TEST	Approved tube types for laboratory testing
<ul style="list-style-type: none">▪ Immunodeficiency Panel (IMD)▪ Lymph Phenotype (TBNK)▪ CD3 Transplant▪ CD19	 <p data-bbox="657 598 738 640">EDTA</p>
<ul style="list-style-type: none">▪ Flow Peripheral Blood	 <p data-bbox="657 1039 738 1081">EDTA</p> <p data-bbox="812 1039 933 1123">Sodium Heparin</p>