

UCD DOPLM TOXICOLOGY ALTERNATE SPECIMEN TUBES

TEST	Approved tube types for laboratory testing
<ul style="list-style-type: none"> - Methylmalonic Acid (MMA) 	<div data-bbox="857 281 932 554"></div> <div data-bbox="857 554 932 596">SST</div> <div data-bbox="1062 281 1127 554"></div> <div data-bbox="1029 554 1192 596">RED TOP</div>
<ul style="list-style-type: none"> - Acid/Neutral Drug Screen - Pentobarbital 	<div data-bbox="857 680 932 953"></div> <div data-bbox="818 953 980 995">RED TOP</div> <div data-bbox="1045 659 1127 974"></div> <div data-bbox="1045 995 1127 1016">EDTA</div>
<ul style="list-style-type: none"> - Ethylene Glycol - Volatile Alcohols 	<div data-bbox="850 1058 932 1373"></div> <div data-bbox="850 1373 932 1394">EDTA</div> <div data-bbox="1045 1100 1127 1331"></div> <div data-bbox="1013 1331 1175 1352">Sodium Fluoride</div>