

UCD DOPLM Tube Options for Transfusion Services

Test	Approved tube types for laboratory testing
<ul style="list-style-type: none"> • Type and Screen • Crossmatch • Phenotype • A1 Lectin • Antibody Titer 	<div data-bbox="771 577 852 976"> </div> <div data-bbox="1039 577 1120 976"> </div> <div data-bbox="771 1018 852 1060">EDTA</div> <div data-bbox="1039 1018 1120 1060">RED TOP</div>
<ul style="list-style-type: none"> • All Ortho Vision Testing (Automated Instrument Testing) • DAT • Fetal Maternal Hemorrhage (FMH) • Eluate 	<div data-bbox="755 1092 836 1491"> </div> <div data-bbox="755 1501 836 1541">EDTA</div>