












Swab Collection Guidelines

C. diff Surveillance Test	MRSA Surveillance Test	Bacterial Culture Swabs -aerobic bacterial culture - fungal culture
<p>Use BLUE BBL CultureSwab only</p> 	<p>Use BD ESwab</p> 	<p>Use BD ESwab</p> 
<p>No RED BD Copan Transystem Swabs</p> 	<p>Temporary alternatives: RED BD Copan Transystem Swabs</p> 	<p>Temporary alternatives: RED BD Copan Transystem Swabs</p> 
<p>No GREEN HealthLink TransPorter Swabs</p> 	<p>GREEN HealthLink TransPorter Swabs</p> 	<p>GREEN HealthLink TransPorter Swabs</p> 
<p>No BD ESwab</p> 	<p>No BLUE BBL CultureSwabs</p> 	<p>No BLUE BBL CultureSwabs</p> 