








## UCD DOPLM MOLECULAR ALTERNATE SPECIMEN TUBES

| TEST   | Approved tube types for laboratory testing  |
|--|---|
| <ul style="list-style-type: none"> <li>● FISH                             <ul style="list-style-type: none"> <li>○ CLL</li> <li>○ Direct harvesting of Peripheral Blood or Bone Marrow</li> </ul> </li> </ul>                  | <div style="text-align: center;">  </div> <p style="text-align: center;">Sodium Heparin      EDTA</p>          |
| <ul style="list-style-type: none"> <li>● Factor V, Factor II and MTHFR mutations</li> <li>● Hereditary Hemochromatosis</li> <li>● Genotyping (CF, HCV)</li> <li>● Viral Loads (HCV, HIV-1, CMV, EBV, HBV, BK-Blood)</li> </ul> | <div style="text-align: center;">  </div> <p style="text-align: center;">EDTA Only</p>                        |
| <ul style="list-style-type: none"> <li>● Human Papillomavirus, High Risk</li> </ul>  | <p style="text-align: center;">ThinPrep PreserveCyt Media Only</p> <div style="text-align: center;">  </div> |
|  | <p style="text-align: center;">See next page</p>  |

## UCD DOPLM MOLECULAR ALTERNATE SPECIMEN TUBES

|   |   |
|---|---|
| <ul style="list-style-type: none"><li>• CT/NG</li></ul>   | <p>COBAS Media Only</p>  |
| <ul style="list-style-type: none"><li>• COVID-2019 RNA Qual</li><li>• RVP</li><li>• Bordetella pertussis</li><li>• HSV (Swab)</li></ul> | <p>UTM Only</p>        |
|   | <p>See next page</p>  |

## UCD DOPLM MOLECULAR ALTERNATE SPECIMEN TUBES

|  |   |
|--|---|
| <ul style="list-style-type: none"><li>• HSV (CSF)</li><li>• Meningitis/Encephalitis Panel</li></ul>  | <p>Sterile Tube Only</p>     |
| <ul style="list-style-type: none"><li>• BK Viral Load - Urine</li></ul>  | <p>Sterile Tube Only</p>  |
| <ul style="list-style-type: none"><li>• Gene Amplification (MDM2, EGFR)</li><li>• Deletion (PTEN, 1p-19q)</li><li>• ALK Rearrangements</li><li>• HER-2</li></ul> | <p>FFPE Only</p>  |