



Sponsored by the Benedictine and Presentation Sisters

Yankton Service Center
 501 Summit Street, Yankton, SD 57078
 Phone: 605-668-8169 OR toll free 866-256-2744
 Fax: 605-668-8168

Request for Supplies

Client: _____ Date: _____
 Address: _____
 Person Requesting: _____ Phone: _____

SUPPLY REQUESTS MUST BE SUBMITTED IN WRITING – BY FAX, COURIER DELIVERY, ETC.

(Expected delivery of supplies within 7 working days.)

Item (unit)	QTY ORDERED	QTY SHIPPED
SPECIMEN COLLECTION:		
Vacutainer-purple K2 EDTA 3 ml (100/flat or ea)		
Vacutainer-blue Citrate NA 1.8 ml (100/flat or ea)		
Vacutainer-blue Citrate NA 4.5 ml (ea)		
Vacutainer-green PST Gel LI Heparin 4.5 ml (100/flat or ea)		
Vacutainer-red NO Gel serum 6 ml (100/flat or ea)		
Vacutainer-gold SST Gel serum 3.5 ml (100/flat or ea)		
Vacutainer-dark green NA Hep 6 ml (ea)		
Vacutainer-Quantiferon TB tubes (ea kit)		
Vacutainer-gray top Na Floride 4 ml (ea)		
Vacutainer-royal blue K2 EDTA 6 ml (ea)		
Vacutainer-royal blue NO additive 6 ml (ea)		
MICROBIOLOGY/VIROLOGY COLLECTION:		
Blood Culture-Aerobic blue (ea)		
Blood Culture-Anaerobic purple (ea)		
Blood Culture-Peds yellow (ea)		
Chlam/GC DNA Amp Probe Collect Kit (female/male) (ea)		
Chlam/GC DNA Amp Probe Collect Kit (Urine kit)(ea)		
eSwab culture swab-aerobic & anaerobic (white) (ea)		
eSwab culture swab-mini tip (green) (ea)		
PCR Influenza and RSV (UTM Media) (ea)		
Media-BAP TSA 5% SB (ea)		
Media-Choc (ea)		
Media-MAC (ea)		
Media-Thio (ea)		
UTM for Ureaplasma/Mycoplasma and virus cultures (ea)		
Ova & Parasite Proto-Fix for O&P with Trichrome (ea)		
Stool Culture-Cairy Blair (C&S) (ea)		
Vaginitis DNA Probe Affirm VP111 (collect kit) (ea)		
innoculation loops 1uL (pk)		
Regan Lowe transport media (ea)		

Item (unit)	QTY ORDERED	QTY SHIPPED
MISCELLANEOUS:		
Requisition forms (ea)		
Fetal Fibronectin (Spec Coll) Kits (ea)		
Glucola 75 or 50 g (ea)		
Labels (Large, Small, or ATLAS type)		
Therapeutic phlebotomy bags		
Diluent 2 (ea bottle)		
Ortho GEL Anti-IgG cards (sleeve)		
Ketone Strips		
Ketone controls		
SPECIMEN TRANSPORT:		
Bio Specimen Transport Bags (6x9)		
Bio Specimen Transport Bags (10x10)		
Transport Vials-Plastic (ea)		
Transport Vials-Plastic (Amber) (ea)		
Transport tube trace element (blue cap)		
URINE COLLECTION:		
24-hr urine container (ea)		
Urine Culture Transport (grey top) (ea)		
UA Tube Transport (Yellow/Red Top) (ea)		
Urine Containers-Sterile Cup (ea)		
CLIENT - or - ASHH COMMENTS/OTHER:		
LAB DOCUMENTATION ONLY		
Order Filled By: _____		
Date Filled: _____		
Sent By: In-town courier Med Express Same Day Other		