


RESPIRATORY INFECTIOUS DISEASES

Winter 2020-2021

For information or questions refer to the online lab test directory, or call the Client Services hotline (24/7): 801-507-2110 or 1-877-353-1106

- 1** Nasopharyngeal flocked swabs in VTM/UTM or
 - 2** Nasal (mid-turbinate) swabs in VTM/UTM or
 - 3** Nasopharyngeal aspirate, suction, or wash
 - 4** Nasal aspirate, suction, or wash
- Throat, mouth, & nose swabs are NOT acceptable



TEST ORDER CODE	TESTS FOR	TAT From RECEIPT in Performing Lab	STABILITY	INTENDED USE	COST System Hospital Labs
FLUPCR^{1,2}	Influenza A Influenza B	RAPID 1-2 hours	Room Temp: 24 hrs Refrigerate: 7 days Frozen: 3 month	ED or other situations where a rapid result will influence treatment	\$\$
FLURSV^{1,2}	Influenza A Influenza B Respiratory Syncytial Virus	RAPID 1-2 hours	Room Temp: 24 hrs Refrigerate: 7 days Frozen: 3 month		\$\$\$
COVQD^{1,2,4}	SARS-CoV-2 Influenza A Influenza B Respiratory Syncytial Virus	RAPID 1-2 hours	Room Temp: 24 hrs Refrigerate: 7 days		\$\$\$
BOPPCR^{1,3}	<i>Bordetella pertussis</i>	STANDARD ≤ next day	Room Temp: 24 hrs Refrigerate: 7 days Frozen: 3 month	Aid to the diagnosis of <i>B. pertussis</i> (whooping cough)	\$\$
RF2PCR¹	Adenovirus Coronavirus 229E Coronavirus HKU1 Coronavirus NL63 Coronavirus OC43 SARS-CoV-2 Metapneumovirus Rhinovirus/Enterovirus Influenza A Influenza B Parainfluenza 1 Parainfluenza 2 Parainfluenza 3 Parainfluenza 4 Respiratory Syncytial Virus <i>B. Pertussis</i> <i>B. parapertussis</i> <i>C. pneumoniae</i> <i>M. pneumoniae</i>	RAPID 2 -3 hours	Room Temp: 4 hrs Refrigerate: 3 days Frozen: 1 month	ED or other situations, including inpatient, where a rapid result (other than Influenza) will influence treatment	\$\$\$\$

Nasopharyngeal (and nasal) swabs are performed at all Intermountain laboratory facilities. All other source types are performed at Central Laboratory or Primary Children's.