

Southern Region Laboratory - DOWNTIME LAB TEST REQUEST

Instructions: Use this form to request laboratory tests when computer systems are down. Inclusion of phone, fax, and tube station numbers is essential for communication of patient test results. In an extended downtime, the Laboratory will limit testing to the items listed below. Please hold requests for other tests until systems are restored. Tests that are not on this list may be written in and will be considered on a case-by-case basis by laboratory directors.

RN Phone Number: _____ Unit Fax Number: _____ Unit Tube Station: _____ Unit/Room # _____		Label with Priority <input type="checkbox"/> STAT <input type="checkbox"/> Timed <input type="checkbox"/> Routine	Patient Information <input type="checkbox"/> Unit Collect: <input type="checkbox"/> Lab Collect
Ordering Physician (Full Name): _____ <input type="checkbox"/> Blood <input type="checkbox"/> Swab <input type="checkbox"/> Urine <input type="checkbox"/> Other (SPECIFY) _____			

***All collection information required**

Collection date: _____	Time: _____	Collected by (full name): _____
-------------------------------	--------------------	--

Order only tests medically necessary for the treatment or diagnosis of the patient. For outpatients, enter the ICD10 codes which demonstrate the medical necessity of each test ordered (REQUIRED)	ICD10 Code(s) _____
--	----------------------------

Chemistry Testing - Panels			Drug Testing			Microbiology/Molecular cont.		
LAB15	Basic Metabolic Panel	GG	LAB3043	Acetaminophen	GG	LAB8619	Blood Culture	BC
LAB17	Complete Metabolic Panel	GG	LAB21	Carbamazepine	GG	LAB268	CSF Aerobic Culture w/ Draw Site(s):	SC
LAB16	Electrolytes Panel	GG	LAB23	Digoxin	G		Gram Stain	
LAB20	Hepatic Function Panel	GG	LAB3051	Ethanol	GG	LAB2022	FLU/RSV/SARS PCR - Nares	VM
Chemistry Testing - Individual Tests			LAB3601	Phenobarbital	G	LAB1895	SARS-CoV-2 PCR	NVM
LAB45	Albumin	GG	LAB31	Phenytoin, Total	GG	LAB8664	Grp A Strep PCR - Throat	ES
LAB3014	Alkaline Phosphatase	GG	LAB3049	Salicylates	GG	LAB1988	Resp Pathogen PCR	NPVM
LAB132	ALT (SGPT)	GG	LAB24	Valproic Acid	GG	LAB2830	Multiplex Vag Panel PCR	PS
LAB4216	Amnisure		LAB3289	Urine Drugs of Abuse	U	Other Culture (list test & body site): Culture/PCR Test: _____ Body Site, Laterality: _____		
LAB131	AST (SGOT)	GG	Coagulation			Transfusion Services		
LAB4413	Bilirubin Direct	GG	LAB320	PT with INR	B	LAB895	ABORH □ Cord Blood (Indicate on tube)	P or SP
LAB50	Bilirubin Total	GG	LAB325	APTT	B	LAB274	DAT □ Cord Blood (Indicate on tube)	P or SP
LAB4311	NT-ProBNP	GG	LAB313	D-Dimer	B	LAB4697	Retype	P
LAB53	Calcium	GG	LAB314	Fibrinogen	B	LAB276	Type and Screen	P
LAB4229	Calcium Ionized	GG	LAB306	Factor VII (STAT only)	B	LAB3494	RhIG Evaluation	P
LAB55	Carbon Dioxide	GG	LAB308	Factor IX (STAT only)	B	LAB4695	Type and Screen Newborn	P
LAB59	Chloride	GG	LAB317	Hep Hybrid (Anti Xa)	B	LAB893	Transfusion Reaction Investigation	P
LAB106	Cardiac BNP (GVH only)	L	LAB316	Low Mol Wt Heparin	B	LAB8370	FMH	P
LAB62	Creatine Kinase (CK)	GG	Hematology and CSF/Body Fluids			*Blood Bank will order test products, copy and give to Blood Bank		
LAB3013	Creatinine	GG	LAB210	CBC w/ Diff	L		Prepare RBC	
LAB82	Glucose Random	GG	LAB294	CBC w/o Diff	L		Specify Qty: _____ □ Irradiated	
LAB85	GGT	GG	LAB322	ESR	L		Prepare Platelets	
LAB3451	HCG Quantitative	GG	LAB287	Fetal Fibronectin (FFN)			Specify Qty: _____ □ Irradiated	
LAB473	HIV 1/2 **See Below	G	LAB8271	CSF Cell Count & Diff	SC		Prepare Plasma, Thawed	
LAB8028	Lactate Plasma	GG	Body Fluid Type:				Specify Qty: _____ □ Irradiated	
LAB8641	Lactate Plasma	GG	LAB8269	Body Fluid Count & Diff	SC		Prepare Cryo, Thawed	
LAB99	Lipase	GG	LAB4949	Crystals, Body Fluid	SC		Specify Qty (4-unit Pools): _____	
LAB103	Magnesium	GG	Urine Evaluation - Specify Source			Other Tests (will be considered by Lab Director)		
LAB4222	Neonatal Bili	GG	□ Clean catch □ Single cath □ Indwelling cath					
LAB113	Phosphorus	GG	LAB4438	UA Complete	U			
LAB114	Potassium	GG	LAB239	Urine Culture	GU			
LAB118	Protein Total	GG	Collection Method:					
LAB122	Sodium	GG	LAB437	Pregnancy Test, Urine	U			
LAB747	Troponin	GG	Microbiology/Molecular					
LAB140	Urea Nitrogen (BUN)	GG	LAB2877	C diff PCR - Stool	SC			
Blood Gases			LAB2005	GI PCR - Stool	SC			
LAB56	Carboxyhemoglobin	GN or S	LAB5200	Mening Enceph Panel	CSF			
LAB76	Arterial Blood Gas	S	LAB4957	GC RNA Source:	AS or SC			
Patient Temperature:			LAB4239	CT RNA Source:	AS or SC			
LAB79	Venous Blood Gas	S						

Unless this testing meets an exception under Colorado law, by authorizing this order you understand that Colorado law requires you to inform the patient that (1) you have ordered testing for sexually transmitted infections; (2) the results may be reported to Colorado's health department; and (3) the patient can opt out of testing.

Memorial Hospital Laboratory - DOWNTIME LAB TEST REQUEST

AMA Approved Chemistry Panels:

Electrolytes Panel	Carbon dioxide, Chloride, Potassium, Sodium
Basic Metabolic Panel	Calcium, Carbon dioxide, Chloride, Creatinine, Glucose, Potassium, Sodium, Urea nitrogen
Comprehensive Metabolic Panel	Albumin, Alkaline phosphatase, ALT (SGPT), AST (SGOT), Urea Nitrogen (BUN), Bilirubin Total, Calcium, Carbon dioxide, Chloride, Creatinine, Glucose, Potassium, Protein Total, Sodium
Hepatic Function Panel	Albumin, Alkaline phosphatase, ALT (SGPT), AST (SGOT), Bilirubin Total, Bilirubin Direct, Protein Total

URINE DRUGS OF ABUSE SCREEN Panel includes Qualitative testing for:

- | | |
|-----------------------|-----------------------------------|
| - Amphetamines | - Methamphetamines |
| - Benzodiazepines | - Opiates |
| - Cocaine Metabolites | - Oxycodone |
| - Fentanyl | - Propoxyphene (PPX) |
| - Marijuana (THC) | - Phencyclidine (PCP) |
| - Methadone | - Tricyclic Antidepressants (TCA) |

Specimen Collection Container Legend:

AS - Hologic Aptima Swab	NVM - Viral Transport Media - Nasal
B - Blue top (3.2% sodium citrate), fill 100%	P - Pink top (EDTA)
BC - Blood culture bottles	PS - Pink capped Xpert Vaginal Swab
CSF - CSF sterile tubes	R - Red top, no gel, clot activator
ES - E-Swab	S - Plastic blood gas syringe with dry balanced heparin
G - Gold top w/ gel (SST clot activator)	SC - Sterile container
GG - Green top w/ gel (lithium heparin)	SP - Short Pink top (EDTA)
GN - Green top, no gel (sodium heparin)	U - plastic urine container
GU - Grey top urine boric acid tube	VP - Pink capped Cepheid Xpert Swab Kit
L - Lavender (purple) top (EDTA)	VM - Viral Transport Media
NP VM - Viral Transport Media - Nasopharyngeal	