



SUPPLY REQUEST FORM

Return Fax Number _____

| |
|------------------------|
| Date: |
| Office: |
| Address: |
| Contact Name/#: |

TriHealth provides collection materials required for the proper collection and transport of specimens that will be submitted to TriHealth's Laboratory Services.
In accordance with the Office of Inspector General's (OIG) Model Compliance Plan for Clinical Laboratories, TriHealth does not provide supplies for the collection of tests that will be referred to other laboratories or performed in the physician's office.

| Quantity | Units | Item |
|--|--------------|--|
| Vacutainer Tubes | | |
| | Box of 100 | Yellow tops (Serum separator tubes) |
| | Box of 100 | Lavender tops |
| | Box of 100 | Blue tops |
| | Box of 100 | Pink tops |
| | Each (glass) | <input type="checkbox"/> Red tops (glass) <input type="checkbox"/> Green tops (lithium heparin) <input type="checkbox"/> Green tops (sodium heparin) |
| | Each | <input type="checkbox"/> Yellow ACD <input type="checkbox"/> Pearl tops <input type="checkbox"/> Tan tops |
| | Each | <input type="checkbox"/> Navy tops (with EDTA) <input type="checkbox"/> Navy tops (without EDTA) |
| Needles | | |
| | Box of 48 | Vacutainers, 21g x 1-1/4" |
| | Box of 48 | Vacutainers, 22g x 1-1/4" |
| | Each | Needle holders (adapters) |
| | Box of 25 | Tourniquets |
| Cytology | | |
| | Tray of 25 | Thin Prep Pap smear vials |
| | Bag Included | Thin Prep Pap spatulas, and cytobrushes |
| | Bags of 25 | Thin Prep Pap brooms |
| Histology | | |
| | Each | 10% formalin containers: <input type="checkbox"/> 15 mL <input type="checkbox"/> 60 mL |
| Microbiology and Virology | | |
| | Set | Blood culture media |
| | Each | Blood culture collection adapters |
| | Each | eSwabs: <input type="checkbox"/> Standard white cap <input type="checkbox"/> Blue cap (for pertussis) |
| | Each | Red capped culturettes (For Grp A Strep Probes only) double flocced |
| | Each | Amplified probe media (GC/Chlamydia): <input type="checkbox"/> Genital <input type="checkbox"/> Urine |
| | Kit | Stool collection: <input type="checkbox"/> 4-vial kit for bacteria, parasites, and/or C. Dife <input type="checkbox"/> 1-vial kit for C. Dife only |
| | Each | Virus transport media (Herpes/Chlamydia) |
| | Each | Quantiferon kits for TB |
| Containers and Kits for Urine and Non-Blood Specimens | | |
| | Each | 4 oz. sterile containers |
| | Box of 50 | Urine culture collection kits (cup, culture tube, and urinalysis tube) |
| | Each | 24-hour urine collection containers with white shopping bags |
| | Each | Nuns cap specimen collection bowl |
| | Box of 100 | Urine transport tubes (yellow top urine tube) |
| | Bag of 100 | Urine transfer straws |
| | Box of 50 | Urine transfer straws with urine culture tubes |
| | Box of 100 | Urinalysis tubes (red and yellow top) |
| Order Forms | | |
| | Pad of 100 | Test requisitions |
| | Pad of 100 | Pap Smear/Biopsy requisitions |
| | Pad of 50 | Supply Request forms |
| Other Items | | |
| | Bottle | Glucose drink: <input type="checkbox"/> Lemon Lime-100 g (only) <input type="checkbox"/> Orange - 50g, 75g, 100g |
| | Pkg of 100 | Specimen transport bags |
| | Pkg of 50 | Large bags |