

COVID Testing: as of 5/20/2020 subject to change

Test Code	Test Name	Source
COVID	SARS CoV2/COVID 19	NP or OP VTM or Saline
COVRT	SARS COV2, COVID 19 RAPID *Available only with preapproval from Dr. Gewirtz or specific ordering locations	NP VTM or Saline

COVID Collection Kits: as of 5/20/2020 subject to change

Collection Kit Type	Supply Location
OP Saline	Materials Warehouse SAP #238127
NP VTM (for use with COVRT Rapid COVID)	IPAC room C-204
OP VTM (For Chandler ED, Drive Thru, Good Sam ED, Turfland Clinic Lab)	Microbiology Lab will distribute upon request

- VTM must be stored refrigerated
- We are no longer performing the Abbott ID NOW 15 minute test
- The following locations have preapproval for COVRT Rapid COVID Testing:
 - Cath lab
 - Endoscopy
 - Transplant
 - OR cases – Trauma or unexpected arrival without previous test
 - Radiation/brachytherapy
 - ECMO
 - Interventional radiology
 - ED transfer to GSH
- If you are concerned about the request for any reason, Dr. Amy Gewirtz can be reached at 614-738-9785.