

Complement fixation (Default

if no order selected)

(+ \$44 if

reflexed)

2 mL

86171

Ship specimens frozen or w icepack (49 CFR 173.199 or IATA PI 602 compliant):

(530) 752-1757 FedEx, UPS, Coccidioidomycosis Serology Lab, Tupper Hall Room 3137, UC Davis, (530) 752-6813 Fax or Courier West Health Sciences Dr., Davis CA 95616 Email coccy@ucdavis.edu Coccidioidomycosis Serology Lab, Tupper Hall Room 3137, UC Davis, USPS www.ucdmc.ucdavis.edu/medmicro/cocci.html 1 Shields Avenue, Davis CA 95616 PATIENT NAME (LAST, FIRST, MI) ORDERING PHYSICIAN (LAST, FIRST, MI) NPI# CLIENT PATIENT ID (MRN) ORDERING PHYSICIAN EMAIL GENDER ICD-10 DIAGNOSIS CODE(S) DATE OF BIRTH F M Other SEND REPORT TO (FACILITY NAME, ADDRESS, STATE, ZIP, EMAIL) ENTITY TO BE BILLED (IF DIFFERENT THAN ORDERING FACILITY) SEND REPORT TO (FAX + OPTIONAL EMAIL OR WEB ACCOUNT) E-FAX COMPATIBLE? YES NO NO SELF-PAY/PATIENT/INSURANCE BILLING SPECIMEN INFORMATION **COLLECTION DATE** SPECIMEN ID# SPECIMEN TYPE Serum Other: Cerebrospinal fluid Lumbar Cisternal 3 mL (MINIMUM 2mL) SPECIMEN REQUIRED WITH REQUESTS FOR IMMUNODIFFUSION AND COMPLEMENT FIXATION 1. Ventricular (VP Shunt) TUBES MUST BE LABELED WITH AT LEAST 2 PATIENT IDENTIFIERS SPECIMENS MUST BE SHIPPED IN COMPLIANCE WITH 49 CFR, PART 173.199 OR IATA PACKING INSTRUCTION 602. Onset of Illness SPECIMEN STABILITY: 6 MONTHS FROZEN, 2 WEEKS AT 4°C Ethnicity TOTAL WBC Lymph_ DIFF: Neut Mono Eos CSF: Glucose_ Protein_ Cells ___ Histo___ Other Immunocompromising conditions/drugs Coccy serologic results obtained elsewhere: ID____Latex_ EIA:(IgG) Coccidioidal skin test____ EIA:(IgM) Brief History/Current patient status: Anti-fungal Medications (include total dose to date and route): **TEST(S) REQUESTED** CPT MIN. DAYS **TEST INFORMATION TEST** COST TAT CODE(S) VOLUME PERFORMED Macrotiter tube-based quantitative IgG assay. In some situations, Coccidioidal complement 2-3 \$44 86171 1 mL quantitative immunodiffusion may also be performed at no M - F fixation (CF) days additional cost. Coccidioidal CF titers correlate with prognosis. Qualitative IgG and IgM. Specimens are concentrated prior to testing, Coccidioidal immunodiffusion 2-3 \$36 86331 1 mL increasing sensitivity. Coccidioidal immunodiffusion is more sensitive M - F days and specific than complement fixation. Reflex: Immunodiffusion to \$36 86331 Immunodiffusion is performed and complement fixation is only 2-5

performed (and billed) if immunodiffusion is positive.

NPI#

CLIA#

Director NPI#

Lab Phone

1386895068

05D0644125

1831206812

days

M - F