

Aspirus Reference Laboratory Client Supply Order Form

Client Name: _____ Date Ordered: _____

Ordered By: _____ Dated Filled: _____

FAX order sheet to 715-847-2576

Note: Supplies ordered must be used only for specimens sent to Aspirus Reference Laboratory

Capillary Collection:

Quantity = **Each**

- _____ Amber (20.1344)
- _____ Mint (20.1345.100)
- _____ Lavender (20.1341.102)
- _____ Red (20.1343.100)
- _____ SST-Amber/Yellow (365978) Bili
- _____ Gray (20.1307.100)

Blood Tubes: Quantity - Indicate tray or each

- _____ Light Blue 1.8ml (363080)
- _____ Blue 2.7 ml (363083)
- _____ Gray 2ml (367921)
- _____ Lavender 4ml (367861)
- _____ Navy Serum 6ml (368380)
- _____ Navy EDTA 6ml (368381)
- _____ Tan 3ml (367855) (lead)
- _____ Green NA Hep 4.0ml (367871)
- _____ Mint Lith Hep.4.5ml (367962)
- _____ Pink 6ml (367899)
- _____ Red 4ml (367812)
- _____ Yellow SST 5.0ml (367986)
- _____ Yellow ACD 8.5ml (364606)

Forms: (Forms can be found in ARL Reference Manual)

Send Outs: Quantity = Each only

- _____ Protease, Refrigerator (49662)
- _____ Saliva Cortisol (52056)
- _____ Sulfamic Acid (5-57300)
- _____ Calculi, Risk/Urines (KT-1080)
- _____ SPS Sterile Tubes (364960)
- _____ Cell-Free DNA Blood Col. Tubes (52358)
- _____ NMR Tube (456073P) - Lipomed
- _____ Natera Test
- _____ TB Gold (ASM02085)
- _____ DIBT Tubes (222-3010-080)
- _____ Serotonin, Asorbic Acid (**only get when DR. orders**)
- _____ Blue Capped Transf Tubes (62.612.014)-Lead Free
- _____ Red Capped Sterile Tubes (62.612.015)-Body Fluid
- _____ 10cc Sterile Tubes (62.551.20)
- _____ Millennium Kit PGT (oral swab kit)
- _____ Millennium Kit UDT (urine kit)
- _____ UTM, (Chlamydia, Myscplasma, Virus Culture) (KT-8800)
- _____ RPMI (1640) (Flow area)
- _____ Body Fluid Kits (IVR/ASP Only)
- _____ Blood Bands/Envelopes (4R4588)

Quantity= **Covid-19 Kits:By The Each**

- _____ NP Swab COVID PCR Collection Kits
- _____ Nasal Swab COVID PCR Collection Kits
- _____ NP Swabs for COVID PCR (**Flexible Mini Tip**)
- _____ Nasal Swabs for COVID PCR (**Rigid Shaft**)
- _____ COVID PCR Transport Media
- _____ Pink COVID Labels
- _____ Yellow COV2AG Labels
- _____ Employee Collection Kit Components (**Symptomatic**)
- _____ **Unvaccinated Employee Screening Components**
- _____ **Aspirus Clinical Laboratory Requisitions (Lab-019)**

Microbiology Supplies:

Quantity = **Box/Each**

- _____ E-Swabs Regular, White Cap (220245)
- _____ E-Swabs Mini, Blue Cap (220532)

Microbiology Supplies:

Quantity = **Box/Each**

- _____ NP Collection/viral transport Swabs (Swab/B100 influenza)
- _____ Blood Culture BD Grey, Aerobic (442023)
- _____ Blood Culture BD Purple, Anaerobic (442021)
- _____ Blood Culture BD Pink, Ped (442020)
- _____ Pertussis PCR NPG (12300/12310) Bactiswab (**Backordered**)
- _____ CT/GC PCR, Aptima Multitest Swab Specimen (PRD-03546)
- _____ CT/GC PCR, Aptima Urine (301040)
- _____ Aptima Specimen Transfer Kit (301154C)
- _____ Para-Pak Clean, White Cap (900312)
- _____ Para-Pak C&S, Orange Cap (900612)
- _____ Para-Pak Formalin, Pink Cap (900412)
- _____ Parasep SF, Green Cap (149995)
- _____ Affirm VP III ATTS (446255)
- _____ Urine, C&S Transport, Grey Cap (364953)
- _____ CO2 Pack, Mitsubishi (R681007)
- _____ CMG Broth (See Micro) (R060544)
- _____ LIM Broth (See Micro)
- _____ Pin Worm Paddles (220120)
- _____ Slide Holders (2 Slides) (423787)
- _____ Anaerobic Tissue Transport (AS-919)
- _____ GasPak EZ Araerobe Gas Generating Pouch (4360683A)

Urine Preservatives: 24 Hr Urine Containers

Quantity = **Each**

- _____ Hydrochloric Acid (HCl) 6N (CSR'S)
- _____ Sodium Hydroxide (NaOH) 5% (CSR'S)
- _____ 24 HR Urine Container (empty)

Specimen Containers:

- _____ Arup Standardized Trans. Tubes (Clear) (100/bag) (62.612.016)
- _____ Arup Standardized Amber Tubes (62.612.017) 10 ea/bag
- _____ Serum Vials With Caps (Clear) 500/Bag (459000)
- _____ **Semen cups pre-weighed (post vasesc/Infertility spec.)**
- _____ Sterile 4oc Containers (Urine, Sputum) (DYND30364)
- _____ Stool Containers (72 Hour Fecal Fat) (T808160-2)
- _____ FOBT Green Cap Vials (OCS1)
- _____ FOBT Mailers (OCPU)

Biohazard Bags: (Aspirus Reference Laboratory)

- _____ **Plain Clear Bags for Making Specimen Collection Kits**
- _____ Blue (Frozen) (RQC 622980)
- _____ Clear (Refrigerate) (RQC 622981)
- _____ Green (Ambient) (RQC 622982)
- _____ Yellow (Pathology) (RQC 622983)
- _____ Red (Stat)
- _____ Specimen Bag /each (12x15) (4812155B)
- _____ Patient Belonging bags

Items Billed to Account:

- _____ Needles: 22g 1.25" (368608)
- _____ Needles: 21g 1.25" (368607)
- _____ Butterfly Needle Set 21g (50/pk) (972112)
- _____ Butterfly Needle Set 23g (50/pk) (972312)
- _____ Butterfly Needle Set w/holder 23g (25/pk) (368653)
- _____ Butterfly Needle Set 25g (25/pk) (972512)
- _____ Female Blood Transf. Device (364880)
- _____ Male Luer Lok Access Device (364902)
- _____ Luer Adapter (100/pk) (367290)
- _____ Media-Blood Agar Plates (10/pk) (R01202)
- _____ Needle Holders (100/bag)or (250/bag) (364815)
- _____ Labels (Roll) - Zebra White (8792)
- _____ Labels (Roll) - Zebra Yellow (26407)
- _____ Labels (Roll) - Micro (07112-0005)
- _____ Ribbons (Roll) - Micro (R-04-00-0041-03)
- _____ Tourniquets (100/pk) (VTK0308)
- _____ Chlora Prep (20/bx) (260299)