

BELLIN HEALTH OUTREACH LABORATORY
SUPPLY WORKSHEET

FACILITY: _____ DATE: _____ FILLED BY: _____

SUPPLY	AMOUNT	UNITS (CIRCLE)
COLLECTION TUBES		
LAVENDER 4.0 ML		BOX (100) EACH
BLUE 1.8 ML		BOX (100) EACH
GREEN PST 3.0 ML		BOX (100) EACH
GOLD SST 4.0 ML		BOX (100) EACH
RED PLAIN 6.0 ML		EACH
PINK BLD BANK EDTA		EACH
MICROTAINER LAVENDER		EACH
MICROTAINER GOLD		EACH
MICROTAINER GREEN		EACH
URINE CULTURE PRESERVATIVE		BOX (50) EACH
URINALYSIS PRESERVATIVE		BOX (50) EACH

SUPPLY	AMOUNT	UNITS (CIRCLE)
COLLECTION SUPPLIES		
BD NEEDLE HOLDER		BAG (250) EACH
BD TRANSFER DEVICE-BLUE TIP		EACH
BD TRANSFER DEVICE-PINK TIP		EACH
NEEDLE 22 X 1.25		BOX (48) EACH
NEEDLE 21 X 1.25		BOX (48) EACH
BLD CULT BOTTLE - ANAEROBIC		EACH
BLD CULT BOTTLE - AEROBIC		EACH
BLD CULT BOTTLE - PEDIATRIC		EACH
**BUTTERFLY 21g		BOX (50) EACH
**BUTTERFLY 23g		BOX (50) EACH
BLOOD CULT PREPS		EACH
QUANTIFERON KITS		EACH

SUPPLY	AMOUNT	UNITS (CIRCLE)
MISCELLANEOUS		
PAP KITS THIN PREP		EACH
FROSTED SLIDES		BOX (50)
SPECIMEN LABELS (SEE NOTE)		ROLL

NOTES: _____

SUPPLY	AMOUNT	UNITS (CIRCLE)
FORMS		
OUTREACH REQUISITION		EACH
HISTOLOGY/TISSUE		PACK (100)
PAP DUPLICATE		PACK (100)
CYTOLOGY NON/GYN		PACK (100)

SUPPLY	AMOUNT	UNITS (CIRCLE)
CONTAINERS		
FORMALIN 1 OZ		BOX (12) EACH
FORMALIN 2 OZ		BOX (12) EACH
URINE CUPS (STERILE)		BAG (100) EACH
SPECIMEN TRANSPORT TUBES		BAG (100) EACH
24 HR URINE CONTAINER		EACH
SPEC TRANSPORT BAGS		PACK (100)
TRANSFER PIPETTES		BOX (500)
ALCOR FIX VIAL (GREEN)		EACH

SUPPLY	AMOUNT	UNITS (CIRCLE)
MICROBIOLOGY		
AEROBIC CULTURETTES (RED)		BAG (50) EACH
APTIMA - ORANGE		BOX (50) EACH
APTIMA - YELLOW		EACH
APTIMA - BLUE		EACH
NASAL SWABS		BOX (100) EACH
NASOPHARYNGEAL SWABS		EACH
STOOL TRANSPORT VIALS (ORANGE)		EACH
PINWORM PADDLES		EACH
VIRAL MEDIA		EACH
SPUTUM BULLETS		EACH

ALL SUPPLIES TO BE USED FOR BELLIN PERFORMED TESTING ONLY