

Legacy Laboratory Services

Legacy LabAlert

10/12/2021

An Important Update from Legacy Laboratory Services

Due to a national reagent backorder the test Trichomonas vaginalis by PCR will not be available for another month. All orders for Trichomonas vaginalis by PCR will be converted to an ARUP send out test upon arrival in the Microbiology department. The test will be cancelled as a test substitution and a Referral order will be placed for the ARUP test (Trichomonas vaginalis by Transcription Mediated Amplification (TMA), see the below table. At this time only urine specimens will be accepted for this test.

Test Name	Acceptable Specimen	Notes	Stability	Turn-around time
Trichomonas vaginalis by Transcription Mediated Amplification (TMA)	Minimum of 2 mLs of Fresh Urine Transport refrigerated	This test is a send out test to ARUP test code 2005506.	Fresh Urine is stable for 24 hours. Once placed in transport tubes – refrigerated: one month Frozen: one year	Performed Monday, Wednesday, Friday. Reported 1-4 days

For additional information, please contact your account representative, client services or consult our website:

Legacy Laboratory Client Services: 503-413-1234, 877-270-5566, www.legacyhealth.org/labservices