

MUSC Laboratory Supply Request (Non-MUSC Clients)

FACILITY NAME: _____ ADDRESS: _____ _____ _____ PHONE NUMBER: _____	SUBMITTER NAME (PRINT FIRST AND LAST NAME): _____ _____ DATE OF REQUEST: _____ DATE OF NEED: _____
--	--

Important Note - These requested supplies are to be used for specimens sent to MUSC Hospital Lab only. Avoid using these supplies for specimens sent to other reference labs. (i.e. LabCorp, MUSCP-PC(CFC) Lab, Roper, etc). Your cooperation will allow MUSC Hospital Lab to continue providing this service to your location.

Directions: Indicate quantity of supplies needed. Send the request via courier OR fax to 843-792-4896. **Supply orders are generally filled within two business days.** If you have any questions or an urgent request, please call Laboratory Client Services at 843-792-0707 for assistance.

Laboratory Requests and Forms	QUANTITY		Blood Collection - Vacutainers	MSCM#	QUANTITY
Cytopathology Requests	_____	_____	Blue (Sodium Citrate) 2.0 mL (50/Tray)	319	_____
Laboratory Supply Requests	_____	_____	Blue (Sodium Citrate) 3.5 mL (50/Tray)	1636	_____
Maternal Screening/Molecular Pathology Requests	_____	_____	Dark Blue - EDTA (Purple Stripe) (100/Tray)	38095	_____
MUSC Laboratory Requests	_____	_____	Dark Blue - No Additive (Red Stripe) (100/Tray)	49258	_____
Surgical Pathology Requests	_____	_____	Gray (Potassium Oxalate) 2.0 mL (100/Tray)	1638	_____
Misc. Collection Supplies	MSCM#	QUANTITY	Green (Gel, Lithium Heparin) 4.0 mL (100/Tray)	190	_____
Blood Culture - Aerobic	2366	_____	Green (No Gel, Sodium Heparin) 4.0 mL (100/Tray)	38202	_____
Blood Culture - Anaerobic	2367	_____	Lavender/Purple (EDTA) 2.0 mL (100/Tray)	1640	_____
Blood Culture - Pediatric, Aerobic	2368	_____	Lavender/Purple (EDTA) 4.0 mL (100/Tray)	1639	_____
Blood Culture Prep Kits	2628	_____	Pink (EDTA) 2.0 mL (100/Tray)	2618	_____
Chlamydia/GC/Trichomonas Aptima Kits (specify):			Plasma Preparation (PPT, White Pearl Top) (100/Tray)	6790	_____
___ Unisex Swab Kit	1152	_____	Red, Plastic (With Clot Activator) 5.0 mL (100/Tray)	1643	_____
___ Urine Kit	1455	_____	Red, Glass (Rubber Top, No Additive) 10.0 mL (100/Tray)	1648	_____
Eswab for Aerobic and Anaerobic - Regular Tip (50/box)	40185	_____	Serum Separator (Gold SST) 4.0 mL (100/Tray)	189	_____
Eswab for Aerobic and Anaerobic - Mini Tip (50/Box)	30	_____	Serum Separator (Gold SST) 5.0 mL (100/Tray)	6857	_____
Swab, Mini Tip Collection Kit, Flexible, Orange (Sinus)	38141	_____	Yellow (Solution A -ACD) 8.5 mL (100/Tray)	186	_____
Stool Container - 24 hour	N/A	_____			
Universal Viral Transport Media (Red Cap) (50/Box)	39656	_____	Blood Collection – Microtainers (600 uL) (50/Bag)	MSCM#	QUANTITY
Urinalysis Preservative Tube (Yellow/Red Marble)	47978	_____	Amber (for Bilirubin)	37866	_____
Urine Container - 24 Hour Plain Jug	453	_____	Green (Lithium Heparin)	140	_____
Urine Culture Preservative Tube (Gray with Straw)	40313	_____	Lavender (EDTA)	1651	_____
Urine Hats for Urine Collection	318	_____	Red (No Gel/No Additive)	966	_____
Urine Tube, No additive (Yellow Plastic)	51103	_____			
Urine Transfer Device (Straw)	49256	_____	Blood Collection - Accessories and Kits		QUANTITY
Microbiology Media	QUANTITY		Specimen Biohazard Bags		_____
CO2 Gas Generator	_____		Labels for Blood Tubes (Name, DOB, Date, Time)		_____
Chocolate II Agar (10/Pack)	_____		Labels for Cerner Barcode Printer		_____
Inhibitory Mold Agar with Gentamycin (10/Pack)	_____		Pour-off Tubes (Plastic, Screw Cap)		_____
Lowenstein-Jensen Slant	_____		Quantiferon Kits		_____
Page's Amoeba Saline	_____		Salivary Cortisol Kits		_____
Thioglycollate Broth	_____				
TSA II 5% Sheep's Blood Agar	_____		Other Supplies (Not Listed Above):		QUANTITY
Cytology and Pathology Supplies	QUANTITY		_____		_____
95% Ethanol ___ Small ___ Large	_____		_____		_____
Cytolyt (20/Flat)	_____		_____		_____
Formalin Vials - 40 mL	_____		_____		_____
Frosted Slides	_____		_____		_____
Michels Media	_____		_____		_____
Pap Smear Accessory – Brushes/Spatulas (25/Pack)	_____		_____		_____
Pap Smear – ThinPrep Vial (25/Flat)	_____		_____		_____
Slide Holders	_____		_____		_____