



Adventist Medical Center
Laboratory Services

LABORATORY SUPPLIES

Please allow up to a week for delivery

Facility / Clinic

Bacteriology Supplies:

- ___ ea C&S Culture Kit (Culturette)
- ___ ea O & P ECOFIX (1 vial)
- ___ ea Enteric Pathogen Transport (1 vial)
- ___ ea Cobas
 - Urine, Female
- ___ ea Herpes Culture Tubes
- ___ ea Urine Culture & U/A Kits (w/Preserv.)
- ___ Box Fetal Fibronectin
- ___ Box Urine Collection - Straw
- ___ ea Flu Culture/ Influenza Pertussis

Containers:

- Cottage Cheese (non sterile)
- ___ ea 4 oz (150mL)
- ___ ea 16 oz (1 pt)
- ___ ea 32 oz (2 pts)
- ___ ea 83 oz (1 qt)
- ___ ea 172 oz (1 gal)
- ___ ea Sterile All Purpose
- ___ ea 24 Hr. Urine (non sterile)
- ___ ea HAT, Urine

Glucose Tolerance Test:

- ___ ea Lemon-Lime 75 gms
- ___ ea Orange 75 gms
- ___ ea Fruit Punch 50 gms
- ___ ea Lemon-Lime 50 gms
- ___ ea Orange 50 gms

Cytology Supplies:

- ___ pk Brooms (Unimar)
- ___ tray Thin Prep
- ___ tray Sure Path
- ___ pk Cytobrush / Spatulas
- ___ ea Cytolyt (for non Gyn collection)

Lab Requisition Slips:

- ___ ea Lab Pick-up forms (Lab pickup specimen)
- ___ pad Lab Supplies
- ___ shts Lab Labels for specimen ID

Other Lab Supplies:

- ___ ea Vacutainer needle holder
- ___ ea Bag double pouch (Haz. label)
- ___ bg Gauze 2x2 (non sterile)
- ___ ea Sharps waste needle container qt.
- ___ ea Sharps waste needle container 2 gal.

Date

- ___ rl Tape 1 inch
- ___ bx Wipes Alcohol
- ___ Box Towelettes
- ___ Box Breathtek

Anatomic Pathology Supplies:

- Tissue Specimin Bottle w/10% Formalin
- ___ ea small ___ ea large

Needles:

- ___ bx BD Eclipse Vacutainer

Tubes:

- ___ ea Blue (Na Citrate) 4.5ml
- ___ ea Dk Blue w/EDTA
- ___ ea Dk Blue w/o EDTA
- ___ ea Gray (Na Fluoride) 5ml
- ___ ea Lavendar (EDTA) 5ml
- ___ ea Light Green
- ___ ea Pink (BB) 6ml
- ___ ea Red (no additive) 4ml
- ___ ea Yellow (SST) 6ml
- ___ FIT Kit 20 env
- ___ FIT Bottles (50 bottles)