



CLIENT SUPPLY ORDER FORM

Please allow 72 hours for delivery

Deliver supplies to:

Clinic Name: _____
 Address: _____
 Phone: _____
 Contact: _____

| Item # | BLOOD COLLECTION TUBES | Quantity | Units | Item # | MICROBIOLOGY SUPPLIES | Quantity | Units |
|---|---|----------|----------|-----------|--|----------|---------|
| 10939 | Tiger (Serum Separator)-8.5 mL | | Tubes(s) | Varies** | Universal Viral Transport Media | | Tube(s) |
| 10928 | Gold (Serum Separator)-5 mL | | Tubes(s) | 10872 | MRSA PCR Collection Swab | | Swab(s) |
| 10913 | Light Green -(Li Heparin Gel) - 4.5 mL | | Tubes(s) | In Lab | Copan eSwab (Aerobes/Anaerobes) | | Each |
| 10775 | Light Green -(Li Heparin Gel) - 3 mL | | Tubes(s) | 21273** | Genexpert Vag/Cervical Collection /Transport Kit | | Each |
| 10914 | Lavender (EDTA)- 4 mL | | Tubes(s) | 21272** | Genexpert Urine Collection/Transport Kit | | Each |
| 10925 | Lavender (EDTA)- 3 mL | | Tubes(s) | 10766 | Blood culture bottle, Aerobic | | Each |
| 10908 | Dark Green (Na Heparin) - 4 mL | | Tubes(s) | 10767 | Blood culture bottle, Anaerobic | | Each |
| 10768 | Light Blue (Na Citrate) - 2.7 mL | | Tubes(s) | 10765 | Blood culture bottle, Pediatric | | Each |
| 10938 | Light Blue (Na Citrate) - 1.8 mL | | Tubes(s) | | | | |
| 13646 | Gray (K oxalate/Na flouride) - 4 mL | | Tubes(s) | Item # | COLLECTION KITS AND SUPPLIES | Quantity | Units |
| 18642 | Red (Plain) - 10 mL | | Tubes(s) | 1216 | Urine/Stool collection hat, 24 oz (100/Cs) | | Each |
| 10922 | Red (Plain) - 6 mL | | Tubes(s) | In Lab | Fetal Fibronectin Collection Swabs | | Each |
| In Lab ^{NS**} | Royal Blue (EDTA), Trace Element Free-7mL | | Tubes(s) | In Lab | Stool Collection Kit | | Each |
| In Lab ^{NS**} | Royal Blue (No additive)- 7 mL | | Tubes(s) | | | | |
| In Lab ^{NS**} | Quantiferon Gold collection kits | | Kits | Item # | URINE | Quantity | Units |
| | | | | 26971 | UA Urinalysis Transport tubes (Red/Yellow) | | Tube(s) |
| | | | | 26970 | Urine C&S Tubes (Gray Top) | | Tube(s) |
| Item # | SPECIMEN TRANSPORT SUPPLIES | Quantity | Units | In Lab** | 24 hour Urine collection containers | | Each |
| In Lab | Sterile Tube (15 mL conical) | | Each | In Lab | Urine Strainer for Kidney Stones | | Each |
| In Lab | Transfer Tube, Light Protected (Amber, 4mL) | | Each | | | | |
| In Lab | Standard transport tubes (4mL) | | Each | Item # | CYTOLOGY / HISTOLOGY | Quantity | Units |
| In Lab | Transfer Tube, Trace Mineral Free (4mL) | | Each | Pathology | Formalin / Biopsy Bottle - 20 mL | | Each |
| 13498** | Specimen Bags - Small 6x9" | | Each | Pathology | Thinprep Vials w/Purple Papette Brooms | | Each |
| 21059** | Specimen Bags - Large 12x15" (Non-stock) | | Each | | | | |
| In Lab** | Slide sleeves | | Each | Item # | ONLINE REQUISITIONS | Quantity | Units |
| In Lab | Patient Belonging Bag for Stool Kits | | Each | *Online | Add on Requisitions | NA | NA |
| | | | | *Online | Cancel order requisitions | NA | NA |
| | | | | *Online | Teton Pathology Test Requisition | NA | NA |
| **Supply limitations in place. Contact Lab Leadership for allocation information. | | | | | | | |
| ^{NS} Non stock item. Delivery times may be increased | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| OTHER ITEMS NOT LISTED: | | | | | | | |
| | | | | | | | |

Online test menu,specimen requirements, and printable requisitions available at: <https://www.testmenu.com/sjmclab>

EMAIL REQUESTS TO: labsec@stjohns.health (preferred)

FAX TO: 307-739-7323

Current form: March 2023

Lab Use Only: Date order received _____ Date processed: _____