









UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

TEST		Approved tube types for laboratory testing
<ul style="list-style-type: none"> BMP CMP LIPID HFP Renal Panel B-HCG Mg, Phosphorus Amylase CRP, HSCRIP 	<ul style="list-style-type: none"> CK GGT Iron Salicylate Ethanol Osmolality Lipase Albumin <p>(includes all subcomponents of panel)</p>	 <p>PST SST RED TOP</p>
<ul style="list-style-type: none"> TROPONIN PROCALCITONIN MYOGLOBIN CK-MB BHB 		 <p>PST</p>
<ul style="list-style-type: none"> Complement C3 Complement C4 Haptoglobin IGG, IGA, IGM, Prealbumin 	<ul style="list-style-type: none"> PSA Estradiol TSH T4 Digoxin Ferritin 	 <p>SST RED TOP</p>
<ul style="list-style-type: none"> Vitamin D 		 <p>PST SST</p>

UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

<ul style="list-style-type: none"> • Lithium • Theophylline 	 <p>RED TOP</p>
<ul style="list-style-type: none"> • Acetaminophen • Carbamazepine • Gentamicin • Methotrexate • Phenobarbital • Phenytoin • Tobramycin • Valproic acid • Vancomycin 	 <p>Sodium Heparin (Dark Green) RED TOP</p> <p>***NOTE: There are 2 Dark Green Pedi tubes. ONLY Sodium Heparin additive is approved for testing, NOT Lithium Heparin.***</p>
<ul style="list-style-type: none"> • Lactate • Glucose tolerance test 	 <p>Sodium Fluoride</p>
<ul style="list-style-type: none"> • BNP • Ammonia • PTH intact • PTH intraoperative 	 <p>EDTA</p>