









UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

TEST		Approved tube types for laboratory testing
<ul style="list-style-type: none"> BMP CMP LIPID HFP Renal Panel B-HCG Mg, Phosphorus Amylase CRP, HSCRp 	<ul style="list-style-type: none"> CK GGT Iron Salicylate Ethanol Osmolality Lipase Albumin Ammonia <p>(includes all subcomponents of panel)</p>	 <p>PST SST RED TOP</p>
<ul style="list-style-type: none"> TROPONIN PROCALCITONIN MYOGLOBIN CK-MB BHB 		 <p>PST</p>
<ul style="list-style-type: none"> Complement C3 Complement C4 Haptoglobin IGG, IGA, IGM, Prealbumin 	<ul style="list-style-type: none"> PSA Estradiol TSH T4 Digoxin Ferritin 	 <p>SST RED TOP</p>
<ul style="list-style-type: none"> PTH intact PTH intraoperative Vitamin D 		 <p>PST SST</p>

UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

<ul style="list-style-type: none">• Lithium• Theophylline	 <p>RED TOP</p>
<ul style="list-style-type: none">• Acetaminophen• Carbamazepine• Gentamicin• Methotrexate• Phenobarbital• Phenytoin• Tobramycin• Valproic acid• Vancomycin	 <p>Sodium Heparin (Dark Green) RED TOP</p>
<ul style="list-style-type: none">• Lactate	 <p>Sodium Fluoride</p>
<ul style="list-style-type: none">• BNP	 <p>EDTA</p>