


UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

TEST	Approved tube types for laboratory testing
- BMP - CMP - LIPID - HFP - Renal Panel - B-HCG - Mg - Phosphorus - Amylase - CRP - HSCRP - CK - GGT - Iron - Salicylate - Ethanol - Osmolality - Lipase - Albumin (includes all subcomponents of panel)	 PST SST Red Top
- Troponin - Procalcitonin - Myoglobin - CK-MB - BHB	 PST
- Complement C3 - Complement C4 - Haptoglobin - IGG - IGA - IGM - Prealbumin - PSA - Estradiol - TSH - T4 - Digoxin - Ferritin	 SST Red Top
Vitamin D	 PST SST
- Lithium - Theophylline	 Red Top
- Acetaminophen - Carbamazepine - Gentamicin - Methotrexate - Phenobarbital - Phenytoin - Tobramycin - Valproic acid - Vancomycin	 Dark Green Sodium-Heparin Red Top ***NOTE: There are 2 Dark Green Pedi tubes. ONLY Dark Green Sodium-Heparin additive is approved for testing, NOT the Dark Green Lithium-Heparin.***
- Lactate - Glucose tolerance test	 Sodium Fluoride

UCD DOPLM CHEMISTRY ALTERNATE SPECIMEN TUBES

TEST	Approved tube types for laboratory testing
- BNP - Ammonia - PTH intact - PTH intraoperative	EDTA 
HCV Screen	SST 
ED HCV (2 tubes drawn)	BOTH SST  AND EDTA 
HIV Ag/Ab Combo Screen	SST  ***Alternatively Lithium-Heparin, Sodium-Heparine, EDTA, ACD, CPD, CP2D, CPDA, Potassium Oxalate, Sodium-Citrate plasma.***